cleadew MPS

ADVANCED CARE SYSTEM

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OphTECS
Europe

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Stand-alone test bacteria (Standard strain)^A

	<i>P.aeruginosa</i>	<i>S.aureus</i>	<i>S.marcescens</i>	<i>C.albicans</i>	<i>F.sofani</i>
Log reduction value (log/mL)	>4.6	>4.6	>4.7	>4.6	>4.3

(Optectics data)

Acanthamoeba^B


Acanthamoeba castellanii (ATCC 50370)

	Trophozoite	Cyst
Log reduction value (log/mL)	>3.2	>2.2

(Optectics data)


Clinically-isolated bacteria

Planktonic^C

(Optectics data)

Biofilm^D

(Optectics data)

¹**MPDS A:** contains polydronium chloride+myristamidepropyl dimethylamine

²**MPDS B:** contains alexidine dihydrochloride+polydronium chloride

³**MPDS C:** contains polyhexamethylene biguanide hydrochloride

^A**Test method:** In accordance with the stand-alone test, $1.0 \cdot 10^8 - 10^9$ cfu/mL of the test strains are inoculated in the disinfectant and allowed to stand for the period as directed. The viable count is measured after that period.

^B**Test method:** $1.0 \cdot 10^5 - 10^6$ cells/mL of trophozoites and $1.0 \cdot 10^5 - 10^6$ cells/mL of cysts are inoculated in the disinfectant and allowed to stand for the period as directed. The viable count is measured after that period.

^C**Test method:** $1.0 \cdot 10^5 - 10^6$ cfu/mL of the test strains are inoculated in each disinfectant and allowed to stand for the prescribed period of each disinfectant. The viable count is measured after that period.

^D**Test method:** $1.0 \cdot 10^7$ cfu/mL of the test strains are inoculated in a plate and allowed to stand for 24 hours to form a biofilm. Then, each disinfectant is added to the plate and allowed to stand for the prescribed period of each disinfectant. The viable count is measured after that period.

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


Wat is Super Moist Dew (SMD) Technology?

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Evaluation of SMD adsorption^E


^E **Test method:** Silicone hydrogel lenses are immersed in fluorescent-labeled sodium hyaluronate or hyaluronate derivatives for 8 hours, and in ISO-PBS overnight. The volume of each adsorbed ingredient is measured after repeating this procedure 30 times, and the adsorption is observed under a confocal laser scanning microscope.

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Wettability Evaluation^F


^{*} p<0.05, Student's t-test

^{*} **4 MPDS D:** contains polyhexamethylene biguanide hydrochloride + polydronium chloride (with sodium hyaluronate)

^F **Test method:** ISO-PBS was dropped on silicone hydrogel lenses that are repeatedly treated with MPDS D or cleadew MPS 13 times, and the contact angle is measured.

Evaluation of tear stability^G

^{**} p<0.01, Student's t-test (cleadew MPS vs MPDS E)

^F **5** Measured with the Wave front analyzer

^G **6 MPDS E:** contains polyhexamethylene biguanide hydrochloride (with sodium hyaluronate)

^G **Test method:** Silicone hydrogel lenses with **MPDS MPS** or **MPDS E** are used for 2 weeks. Non-invasive break-up time (NIBUT) on the lenses is measured with DR-1 on the morning of day 1, and on the evening of day 3, day 7 and day 14.

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

Clinical studies^H

Efficacy against microbial contamination

Safety for eyes and lenses

^H7 MPDSC C: contains polyhexamethylene biguanide hydrochloride

^H8 ISO 11981, FDA 510 (k)

^HH Test method: With **cleadew MPS**, the subjects wear Group I or IV soft contact lenses for 6 months and 3 months, respectively.

The efficacy (findings of ocular infection and microbiological test) and safety (changes in lens properties and presence of adverse reactions) are assessed.


(Number of cases)

Group I: 172 eyes of 86 patients (safety evaluation), 152 eyes of 76 patients (efficacy evaluation)

Group IV: 80 eyes of 40 patients (safety evaluation), 72 eyes of 36 patients (efficacy evaluation)

^HI Test method: After immersion of the silicone hydrogel lenses in MPDSC C or **cleadew MPS**, the subjects wear them. Two hours later, the condition of the eye surface and total score^I are evaluated. ^IArea x Density Score

Corneal staining test^I

Compatibility of **cleadew MPS** with contact lenses

According to the notification^H, physical, chemical, and biological evaluation were conducted on the lenses after treating Group I and IV lenses 30 times.



The findings showed that **cleadew MPS** can be used for all types of soft contact lenses.

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