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User instructions Microlens GP contact lenses

microlens
you will see

OPTHECS Europe B.V.

Marithaime 20,

6662 WD Elst, The Netherlands

Lens care

- Always wash your hands before touching the contact lens.
- Do not let the contact lens come into contact with soap, lotion, deodorant or any other body product because this can cause irritation of the eye when inserting the lens.
- Never use sharp objects or long fingernails to handle the soft contact lenses.
- Avoid using aerosol products like hairspray or paint spray during lens wearing. If sprays are used wear safety goggles or close the eyes.
- Never use tap water. Always use daily refreshed lens solutions, **never use tap water**.

Description

The GP contact lenses from Microlens are manufactured with ultra-precision computer-controlled lathe cut technology according to the recipe of your contact lens specialist / optometrist / ophthalmologist.

Material

PMMA-silicone co-polymer, Fluorosilicone acrylate, Silicone-hydrogel

Packaging

Microlens GP contact lenses are delivered dry or in a storage solution to your contact lens specialist / optometrist / ophthalmologist. Before insertion the lenses must be kept in a GP contact lens storage solution for at least two hours. This is to optimize the surface properties.

Indications

Microlens GP contact lenses are intended for the correction of visual acuity defects or the optimization of the vision.

Contra-indications

- Inflammations of the anterior eye segment
- Eye diseases that affect the cornea or conjunctiva.
- Insufficient tear secretion
- Systemic diseases that affect the eye.
- When GP contact lenses cause or increase eye allergies.
- Medication with a contra-indication for GP contact lenses (See the instructions for use before taking the medication)

Fitting and wearing instruction

The Microlens GP contact lenses are only fitted, checked and delivered by qualified optometrist / contact lens specialist / ophthalmologist.

Aftercare

- The eyes and contact lens fit should be regularly checked according to the instructions of your optometrist / contact lens specialist / ophthalmologist.
- Always follow the user instructions and prescribed contact lens solutions of your optometrist / contact lens specialist / ophthalmologist.
- Inform your optometrist / contact lens specialist / ophthalmologist about the use of medication which may affect successful contact lens wear.
- Visible damaged lenses should not be worn.
- Always follow the prescribed maximum wearing time.

Possible side effects

In the case of the following sensations:

- burning, stinging or itching eyes;
 - reduced wear time;
 - increased sensitivity to light;
 - excessive tearing, eye discharge or abnormal redness of the eyes;
 - reduced visual acuity;
- Immediately contact your optometrist / contact lens specialist / ophthalmologist.

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